

Multi-Unit Abutment System – RESTORE®

KEY BENEFITS

- Straight and angled (17° & 30°) abutments allow for consistent prosthetics
- Pre-mounted positioning tool for angled abutments makes seating simple
- Designed to accommodate fully and partially edentulous arches
- Torque recommendation: 30 Ncm
- Requires a Quad Driver for all abutments and screws. Recommended torque for prosthetic screws is 20 Ncm



Prosthetic versatility for the multiple unit screw-retained restoration



External Hex connection



External Hex connection



External Hex connection

Implant diameter Ø3.3 mm:

PLATFORM	DESCRIPTION	CUFF	REFERENCE
3.3 mm	straight	1 mm	70413-1
3.3 mm	straight	2 mm	70413-2
3.3 mm	straight	3 mm	70413-3
3.3 mm	angled 17°	3 mm	KM-70408C

Implant diameter Ø3.75 / Ø4.0 mm:

PLATFORM	DESCRIPTION	CUFF	REFERENCE
4.1 mm	straight	1 mm	70311-1
4.1 mm	straight	2 mm	70311-2
4.1 mm	straight	3 mm	70311-3
4.1 mm	angled 17°	3 mm	KM-70238C
4.1 mm	angled 30°	4 mm	KM-70239C

Implant diameter Ø5.0 / Ø6.0 mm:

PLATFORM	DESCRIPTION	CUFF	REFERENCE
5.0 mm	straight	1 mm	70409-1
5.0 mm	straight	2 mm	70409-2

PROSTHETIC COMPONENTS



Description	Reference
Impression Coping, open tray	KM-80016C
Impression Coping, closed tray	KM-80018C
Abutment Analog	KM-90016
Healing Cap	KM-70274
Temporary Sleeve, titanium grade 5, screw included	KM-70157C
Plastic Sleeve, screw included	KM-60030C
Prosthetic Screw	KM-40098
Pre-mounted Positioning Tool	Included with the angled abutments

Items illustrated are not to scale.

*Products may not be available in all markets.